

Jenx Bee Seat

Fully adjustable, modular seating system specifically designed for young children. Versatile range of base, head, trunk, leg and foot support options to precisely tailor the seat to the child's needs.

Playful & gentle introduction to Supported Seating

BRIGHT AND FRIENDLY
thoughtful design for parent and child acceptance


OPTIMUM SEATING POSITION in correct hip angle due to back angle adjustability



ANGLE ADJUSTABLE FOOTPLATE to accommodate or encourage dorsi or plantar flexion

Stable bases provide **GREATER INDEPENDENCE AND INTERACTION** with friends and family

Base Options



Y-BASE: Most flexible base with 25.8-75 cm seat height adjustment. Tilt-in-space of 30° recline - 10° prone.
#3615-8740-000

WOODEN HIGH BASE: Fixed seat height (55 cm). Compact and economical with 20° of Tilt-in-space.
#3615-8720-000

WOODEN LOW BASE: Fixed seat height (40 cm). Compact and economical with 20° of Tilt-in-space.
#3615-8710-000

Tray Options

- Armrests (essential for tray) #3430-5501-000
- Nursery Tray #3430-7101-000
- Tray Elbow Blocks #3430-7350-000
- Activity Frame #3615-7475-000
- Grab Bar #3615-7200-000
- Tray Cut-Out #3430-7111-000



Accessories

Multigrip Headrest
1 Link #3430-4601-000
2 Link #3430-4602-000
3 Link #3430-4603-000

Slimline Laterals
Short #3430-3151-000
Long (with strap) #3430-3251-000




Waistcoat Harness
#3615-6701-000

Shoulder Protractor
#3430-2750-000



Pommel
Size 0 #3615-1650-000
Size 1 #3615-1651-000

Long Leg Sitting Board #3430-0281-000



Knee Blocks
#3430-1751-000

Padded Footplate
#3615-5301-000



Booties
#3430-6065-000

Sandals
#5356-6081-000



CODE	AGE	MAX USER WEIGHT	BACK ANGLE	BACK PAD HEIGHT	CHEST WIDTH	HIP WIDTH	SEAT DEPTH	LOWER LEG LENGTH	DEPTH WITH LONG LEG BOARD
3430-0001-000	4 mths - 5 yrs	25 kg	10° recline - 15° prone	27-36 cm (seat to top)	15-27 cm	13-25 cm	11-35 cm	9.5-34.5 cm	37-56 cm

BEE