

Jenx Corner Seat

Bright, comfortable and ideal first seat for any child, which also suits infants and young children. Offers support in long sitting, which is children's usual first sitting position.

Adjustable floor based seating system

Sturdy wings and corner shape helps to **maintain a midline position** with symmetry of the child's spine

Adjustable Height Backrest can be used with or without the upper back support to support head or shoulders



Abduction Block helps keep the child well positioned, with their hips neatly back into the seat

Small 'touch and close' tape straps keep the Corner Seat **safely connected** to the nursery table while in use

Sitting Benefits

- Long sitting, with the legs stretched out in front, helps to maintain the length of the muscles behind the knee.
- Great working space for peer group interaction and plenty of fun activities including working and playing with other children.



NURSERY TABLE

- Effective when used together with the Corner Seat.
- Large plastic play/activity surface which can be used flat, angled or as an easel.
- Continuous height adjustment to allow the precise and correct height.
- A short strap attaches to the Corner Seat to prevent a child tipping the table over.
- Easy to wipe clean, with cut out to allow the child to get close to the tabletop activities.
- Bright and contrasting edging around the surface reduces the risk of objects sliding off and provides a place for toys to be secured to.



PACKAGE OPTIONS

Package 1 3140-0141-XXX
Corner Seat Size 1 with Nursery Table Size 0

Package 2 3140-0142-XXX
Corner Seat Size 2 with Nursery Table Size 1

Package 3 3140-0143-XXX
Corner Seat Size 3 with Nursery Table Size 2

ACCESSORY

Triangular Cushions

For when a flat back with wings is preferred. Fit all seat sizes. Can also be used beside legs to prevent external rotation.

Lime
3140-2450-054

Blue
3140-2450-066



Size	Code	Approx Age Range	Back Pad Height	Seat Depth	Max User Weight	Base Size	Colour
1	3140-0041-054	6 mths - 2 yrs	31.5 - 42.5 cm	35 cm	25 kg	49 x 50 cm	Lime
2	3140-0042-054	1-6 yrs	44 - 62 cm	37 cm	35 kg	52 x 54 cm	Lime
3	3140-0044-066	4-10 yrs	50 - 68 cm	46 cm	55 kg	60 x 63 cm	Blue