

Jenx Junior +

Comfortable, contoured and stable modular seating system, designed to suit the needs of children with either complex or mild needs. Outstanding levels of adjustment with a wide variety of support options.

Functionality, flexibility and style

Carefully designed to draw attention to the child rather than the equipment itself



Pelvic control system with independently adjustable contoured seating

Flip-up angle adjustable footplate to accommodate dorsi or plantar flexion

All support pads and accessories are designed to quickly & easily flip away without being removed



Y-Base

- Versatile and flexible, with lockable swivel castors
 - High-low range of 25.8-75 cm seat height from floor
 - Tilt-in-space from 30° recline - 10° prone
 - Base footprint of 65 x 78 cm
- 3615-8740-000



Low seat-to-floor height for peer interaction.



Ease of transfer with tilt-in-space and swing-away supports.



Open design for comfort and support.

ACCESSORIES

Multigrip Head Support

- 1 Link 5356-4601-000
- 2 Link 5356-4602-000
- 3 Link 5356-4603-000



Waistcoat Harness

- Size 1 3615-6701-000
- Size 2 3615-6702-000
- Size 3 3615-6703-000



Abduction Block

- 1 3615-1650-000
- 2 3615-1651-000
- 3 3615-1652-000



Shoulder Protractors¹

- 1 3612-2750-000
- 2 3612-2751-000
- 3 3612-2752-000



Lateral Supports¹

- Short 1 3610-3151-000
- Short 2 3615-3152-000
- Flexi 1 3615-3251-000
- Flexi 2 3615-3252-000
- Flexi 3 3615-3253-000




Tray²

- Junior 1 3615-7101-000
- Junior 2 3615-7102-000
- Multiset 3614-7103-000



Tray Infill

- 1 3610-7180-000
- 2 3610-7182-000



Tray Elbow Blocks

- 3615-7350-000



Angle Adjustable Push Handle

- 3612-8610-000



Activity Frame

- 3615-7475-000



Size	Code	Approx Age Range	Max User Weight	Back Angle	Seat Depth	Hip Width	Back Support Height	Back Support Height With Extension	Lower Leg Length	Base Option
1	3612-0001-000	18 mths-5 yrs	50 kg	15° recline - 15° prone	21-30.5 cm	17-37.5 cm	31-41 cm	31-48.5 cm	19.5-27 cm	Y-Base
2	3612-0002-000	3-8 yrs			24-32 cm		33-43 cm	33-50.5 cm	25-36.5 cm	
3	3612-0003-000	5-12 yrs			27.5-39 cm		39-52 cm	39-59.5 cm	25-40.5 cm	
4*	3614-0014-011	8-14 yrs	100kg	-10° - +25°	28.5-45 cm	20-44 cm	46.5-59 cm	N/A	25-47 cm	X-Base
5*	3614-0025-011	10-Adult			35.5-52 cm		52.5-66 cm	N/A		



* Jenx Multiset
1 Requires covers
2 Tray padding options available