

Jenx Prone Stander



**UPRIGHT/PRONE
STANDING SYSTEM**

Client Details

Client/User Name:		Notes
Date:		
NDIS/Funding:		

Size Options

Measurements	Size 1 <input type="radio"/>	Size 2 <input type="radio"/>
Code, \$	5080-0011-000	5080-0021-000
Age Range	3-7 years	6-14 years
Max User Weight	55kg	70kg
Base Size	90 x 65cm	65 x 99.5cm
Chest Width	19-29cm	23-35cm
Hip Width	19-29cm	23-35cm
Hip Pad Height	41-75cm	48-96cm
Chest Height from Footplate	66-100cm	82-127cm
Tilt-in-Space	0-40°	0-40°
Sandals	170mm (L) x 80mm (W)	220mm (L) x 95mm (W)

Package includes:

- Angle-adjustable frame
- Hip and chest pads
- Knee blocks
- Sandals
- Tray with bowl



Accessories

Chest Extension Wedge	Tray Elbow Blocks <input type="radio"/>
<p>Fills gaps in chest area if a user is too small for the frame</p> <p>Size 1 <input type="radio"/> 5080-2931-000</p> <p>Size 2 <input type="radio"/> 5080-2932-000</p>	<p>Encourages arms into a midline position and block excessive retraction</p> <p>5356-7351-000</p>

Sandal Raising Blocks
<p>Gives the sandals that extra bit of height where needed (kit of 4 blocks)</p> <p>Size 1 <input type="radio"/> 5356-6091-000</p> <p>Size 2 <input type="radio"/> 5356-6092-000</p>

medifab®

Distributed by Medifab Ltd

32 Detroit Drive
Rolleston 7675
New Zealand

+64 3 307 9790

solutions@medifab.com

www.medifab.com



DISTRIBUTED BY:

Script: Jenx_Prone_20181026_Unpriced