

Jenx Standz - Prone Option 2

Date: _____ Order Reference Number: _____

Childs name: _____ Child's D.O.B. _____

Therapists Name: _____ Therapists Email: _____

Therapists Phone Number: _____

Therapists Address: _____

Delivery Address: _____

Delivery Phone Number: _____

MEASUREMENTS	APPROXIMATE AGE RANGE 1 - 9 years	Client measurements	Product measurements
	MAX USER WEIGHT		45 KG
	HIP/TRUNK WIDTH		18 - 30 CM
	CHEST HEIGHT FROM FOOTPLATE		58 - 100 CM
	MAXIMUM USER HEIGHT	N/A	138 CM
	INDIVIDUAL LEG ABDUCTION	N/A	0 - 30 degrees
	FRAME ANGLE ADJUSTMENT	N/A	0 - 90 degrees
	BASE SIZE	N/A	64 x 89 CM

STANDZ - PRONE OPTION 2



Standz - Prone Configuration Option 2

Comes with:

- Gas strut assisted angle adjustable base with angle gauge and incorporated safety lock. 4 locking casters.
- Wipe clean yellow vinyl cover set with inbuilt antimicrobial protection.
- 0 - 30 degrees of individual leg abduction with angle gauge.
- Angle adjustable footplates and size 2 sandals.
- Width, height and depth adjustable Knee supports
- Width adjustable pelvic support pads with incorporated de-rotational strap and safety buckle.
- Aysymmetrically adjustable Trunk support with width and height adjustment and incorporated de-rotational strap with safety buckle.
- Tray support bracket and angle adjustable tray with extension wings.
- Pelvic support board.
- Pelvic support board back pad.
- Tool Free Adjustments.

STZ01-P2-09	