

Neat Zeat

Modular indoor seating solution for activity, feeding and learning opportunities. Constructed of durable hardwood with a lacquer finish, the Neat Zeat looks great in any home or school environment.

Family time, Feeding and Play



Extensive growth and adjustability: from **6 months to adolescent**

Adjustable lateral support upgrade available

Angle adjustable seat; choose posterior or encourage an **active sitting posture** in anterior tilt

Huge growth adjustment range which can be easily adjusted using an Allen key

New Zealand made!

Package Options:

Package 1 3320-0001-000

Seat and footplate, padded seat cushion, padded back support, hip belt

Package 2 3320-0002-000

Seat and footplate, padded seat cushion, padded back support, hip belt, adjustable laterals, chest strap

Package 3 3320-0003-000

Seat and footplate, padded seat cushion, padded back support, hip belt, adjustable laterals, chest strap, depth adjustable tray table, mobile base



Angle adjustable seat



Adjustable seat height



Angle adjustable backrest

ACCESSORIES

Mobile Base

3320-8703-000



Tray Table

3320-7102-000




Transfer Step

3320-5251-000



Transfer Handles

3320-5502-000



Laterals

3320-3001-000



Chest Strap

3320-6803-000



Spex Vest Harness

XSm: 1409-6610-017

Sml: 1409-6611-017



Shoulder Protectors

3320-2750-000



Head Support

3320-4611-000



Foot Supports

3320-6075-000



Quick-Release Pommel

3320-1511-000



High Back Support Upgrade

3320-2239-UPG



Code	Approx Age Range	Back Pad Height	Seat Width	Back Height	Footplate Height	Seat Height From Floor	Seat Angle	Back Angle
3320-0001-000	6 months - Adolescent	Kindy: 9-31 cm Std.: 18-36 cm	40 cm	Std.: 28-40 cm High: 40-50 cm	8-56 cm	40 - 68 cm	0°, 5°, 20°	61-115°